
mh-utils

Release 0.2.2

**Utilities for handing ancillary files produced by
MassHunter.**

Dominic Davis-Foster

Apr 25, 2024

Contents

1 Installation	3
1.1 from PyPI	3
1.2 from Anaconda	3
1.3 from GitHub	3
2 mh_utils.utils	5
2.1 as_path	5
2.2 camel_to_snake	5
2.3 element_to_bool	5
2.4 make_timedelta	5
2.5 strip_string	6
3 mh_utils.xml	7
3.1 XMLFileMixin	7
3.2 get_validated_tree	7
4 mh_utils.cef_parser	9
4.1 Compound	10
4.2 CompoundList	12
4.3 Device	13
4.4 Flag	14
4.5 LocationDict	14
4.6 Molecule	15
4.7 parse_cef	15
4.8 parse_compound_scores	15
4.9 parse_match_scores	15
4.10 Peak	16
4.11 RTRange	17
4.12 Score	18
4.13 Spectrum	18
5 mh_utils.csv_parser	21
5.1 ResultParser	21
5.2 parse_masshunter_csv	21
5.3 mh_utils.csv_parser.classes	22
5.4 mh_utils.csv_parser.utils	30
6 mh_utils.worklist_parser	31
6.1 read_worklist	31
6.2 mh_utils.worklist_parser.classes	31
6.3 mh_utils.worklist_parser.columns	41
6.4 mh_utils.worklist_parser.enums	43

6.5	mh_utils.worklist_parser.parser	43
6.6	Example Usage	44
Python Module Index		47
Index		49

The current utilities are as follows:

- *mh_utils.cef_parser*: Parse Agilent MassHunter Compound Exchange Format files (*.cef files).
- *mh_utils.csv_parser*: Parser for CSV result files produced by MassHunter Qualitative.
- *mh_utils.worklist_parser*: Parse Agilent MassHunter Worklists (*.wkl files).

Installation

1.1 from PyPI

```
$ python3 -m pip install mh_utils --user
```

1.2 from Anaconda

First add the required channels

```
$ conda config --add channels https://conda.anaconda.org/conda-forge
$ conda config --add channels https://conda.anaconda.org/domdfcoding
```

Then install

```
$ conda install mh_utils
```

1.3 from GitHub

```
$ python3 -m pip install git+https://github.com/PyMassSpec/mh_utils@master --user
```


`mh_utils.utils`

General utility functions.

Functions:

<code>as_path(val)</code>	Returns <code>val</code> as a <code>PureWindowsPath</code> , or <code>None</code> if the value is empty/ <code>None/False</code> .
<code>camel_to_snake(name)</code>	Convert name from CamelCase to snake_case.
<code>element_to_bool(val)</code>	Returns the boolean representation of <code>val</code> .
<code>make_timedelta(minutes)</code>	Construct a timedelta from a value in minutes.
<code>strip_string(val)</code>	Returns <code>val</code> as a string, without any leading or trailing whitespace.

`as_path(val)`

Returns `val` as a `PureWindowsPath`, or `None` if the value is empty/`None/False`.

Parameters `val` (`Any`) – The value to convert to a path

Return type `Optional[PureWindowsPath]`

`camel_to_snake(name)`

Convert name from CamelCase to snake_case.

Parameters `name` (`str`) – The CamelCase string to convert to snake_case.

Return type `str`

`element_to_bool(val)`

Returns the boolean representation of `val`.

Values of `-1` are counted as `True` for the purposes of this function.

`True` values are `'y'`, `'yes'`, `'t'`, `'true'`, `'on'`, `'1'`, `1`, `-1`, and `'-1'`.

`False` values are `'n'`, `'no'`, `'f'`, `'false'`, `'off'`, `'0'`, and `0`.

Raises `ValueError` if ‘`val`’ is anything else.

Return type `bool`

`make_timedelta(minutes)`

Construct a timedelta from a value in minutes.

Parameters `minutes` (`Union[float, timedelta]`)

Return type `timedelta`

Changed in version 0.1.0: Moved from `mh_utils.cef_parser`.

strip_string(*val*)

Returns *val* as a string, without any leading or trailing whitespace.

Parameters **val** (`str`)

Return type `str`

mh_utils.xml

Functions and classes for handling XML files.

Classes:

<code>XMLFileMixin()</code>	ABC mixin to provide a function for instantiating the class from an XML file.
-----------------------------	---

Functions:

<code>get_validated_tree(xml_file[, schema_file])</code>	Returns a validated lxml objectify from the given XML file, validated against the schema file.
--	--

class XMLFileMixin

Bases: `ABC`

ABC mixin to provide a function for instantiating the class from an XML file.

Methods:

<code>from_xml(element)</code>	Construct an object from an XML element.
<code>from_xml_file(filename)</code>	Generate an instance of this class by parsing an from an XML file.

abstract classmethod from_xml(element)

Construct an object from an XML element.

classmethod from_xml_file(filename)

Generate an instance of this class by parsing an from an XML file.

Parameters `filename` (`Union[str, Path, PathLike]`) – The filename of the XML file.

get_validated_tree(xml_file, schema_file=None)

Returns a validated lxml objectify from the given XML file, validated against the schema file.

Parameters

- `xml_file` (`Union[str, Path, PathLike]`) – The XML file to validate.
- `schema_file` (`Union[str, Path, PathLike, None]`) – The schema file to validate against. Default `None`.

Return type `_ElementTree`

Returns An lxml ElementTree object. When `.getroot()` us called on the tree the root will be an instance of `lxml.objectify.ObjectifiedElement`.

Chapter
FOUR

`mh_utils.cef_parser`

Parser for MassHunter Compound Exchange Format .cef files.

A CEF file represents a file identified in LC-MS data by MassHunter Qualitative. It consists of a list of compounds encapsulated in a `CompoundList`.

A `CompoundList` consists of `Compound` objects representing the individual compounds identified in the data. Each `Compound` object contains information on the location of that compound within the LC data (`location`), the scores indicating the confidence of the match (`compound_scores`), a list of possible matching compounds (`results`), and the matching mass spectrum extracted from the LC-MS data (`spectra`).

The following diagram represents this structure:

- `CompoundList`
 - `Compound`
 - * `Compound.algo` ⇒ str
 - * `Compound.location` ⇒ Optional [`LocationDict`]
 - * `Compound.compound_scores` ⇒ Optional [`Dict [str, Score]`]
 - * `Compound.results` ⇒ List
 - `Molecule`
 - Another `Molecule`
 - ...
 - * `Compound.spectra` ⇒ List
 - `Spectrum`
 - Another `Spectrum`
 - ...
 - Another `Compound`
 - ...

Classes:

<i>Compound</i> ([algo, location, compound_scores, ...])	Represents a compound identified in mass spectral data by MassHunter Qualitative.
<i>CompoundList</i> ([instrument, compounds])	A list of Compound objects parsed from a CEF file.
<i>Device</i> (device_type, number)	Represents the device that acquired a <i>Spectrum</i> .
<i>Flag</i> (string, severity)	Represents a flag in a score, to warn that the identification of a compound is poor.
<i>LocationDict</i>	<i>TypedDict</i> representing the location of a spectrum within mass spectrometry data.
<i>Molecule</i> (name[, formula, matches])	Represents a molecule in a CEF file.
<i>Peak</i> (x, rx, y[, charge, label])	A peak in a Mass Spectrum.
<i>RTRange</i> ([start, end])	Represents an <RTRange> element from a CEF file.
<i>Score</i> (score[, flag_string, flag_severity])	A score indicating how well the compound matches the observed spectrum.
<i>Spectrum</i> ([spectrum_type, algorithm, ...])	Agilent CEF Spectrum.

Functions:

<i>make_timedelta</i> (minutes)	Construct a timedelta from a value in minutes.
<i>parse_cef</i> (filename)	Construct an <i>CompoundList</i> object from the given .cef file.
<i>parse_compound_scores</i> (element)	Parse a <CompoundScores> element into a mapping of algorithms to scores.
<i>parse_match_scores</i> (element)	Parse a <MatchScores> element into a mapping of algorithms to scores.

class Compound(algo='', location=None, compound_scores=None, results=None, spectra=None)

Bases: *Dictable*

Represents a compound identified in mass spectral data by MassHunter Qualitative.

Parameters

- **algo** (*str*) – The algorithm used to identify the compound. Default ''.
- **location** (*Optional[LocationDict]*) – A dictionary of information to locate the compound in the spectral data. Default *None*.
- **compound_scores** (*Optional[Dict[str, Score]]*) – A dictionary of compound scores. Default *None*.
- **results** (*Optional[Sequence[Molecule]]*) – A list of molecules that match the spectrum. Default *None*.
- **spectra** (*Optional[Sequence[Spectrum]]*) – A list of spectra for the compound. Default *None*.

Methods:

<i>__repr__</i> ()	Returns a string representation of the <i>Compound</i> .
<i>__str__</i> ()	Returns the <i>Compound</i> as a string.
<i>from_xml</i> (element)	Construct a <i>Compound</i> object from an XML element.

Attributes:

<code>algo</code>	The algorithm used to identify the compound.
<code>compound_scores</code>	A dictionary of compound scores.
<code>location</code>	A dictionary of information to locate the compound in the spectral data.
<code>results</code>	A list of molecules that match the spectrum.
<code>spectra</code>	A list of spectra for the compound.

`__repr__()`

Returns a string representation of the `Compound`.

Return type `str`

`__str__()`

Returns the `Compound` as a string.

Return type `str`

`algo`

Type: `str`

The algorithm used to identify the compound.

`compound_scores`

Type: `Dict[str, Score]`

A dictionary of compound scores.

`classmethod from_xml(element)`

Construct a `Compound` object from an XML element.

Parameters `element` (`ObjectifiedElement`) – a Compound XML element from a CEF file.

Return type `Compound`

`location`

Type: `LocationDict`

A dictionary of information to locate the compound in the spectral data.

`results`

Type: `List[Molecule]`

A list of molecules that match the spectrum.

`spectra`

Type: `List[Spectrum]`

A list of spectra for the compound.

class CompoundList (instrument='', compounds=None)

Bases: NamedList

A list of Compound objects parsed from a CEF file.

The full `list` API is available for this class.

Parameters

- **instrument** (`str`) – String identifying the instrument that acquired the data. Default ''.
- **compounds** (`Optional[Iterable[Compound]]`) – List of compounds identified in the mass spectrometry data. Default `None`.

Methods:

<code>__repr__()</code>	Return a string representation of the <code>NamedList</code> .
<code>__str__()</code>	Returns the list as a string.
<code>from_xml(element)</code>	Construct a <code>CompoundList</code> object from an XML element.

Attributes:

<code>instrument</code>	The type of instrument that obtained the data, e.g.
-------------------------	---

`__repr__()`

Return a string representation of the `NamedList`.

Return type `str`

`__str__()`

Returns the list as a string.

Return type `str`

`classmethod from_xml(element)`

Construct a `CompoundList` object from an XML element.

Parameters `element` (`ObjectifiedElement`) – The XML element to parse the data from.

Return type `CompoundList`

`instrument`

Type: `str`

The type of instrument that obtained the data, e.g. "LCQTOF".

class Device(device_type, number)

Bases: `object`

Represents the device that acquired a `Spectrum`.

Parameters

- `device_type` (`str`) – String identifying the type of device.
- `number` (`int`)

Attributes:

<code>device_type</code>	String identifying the type of device.
<code>number</code>	

Methods:

<code>from_dict(d)</code>	Construct an instance of <code>Device</code> from a dictionary.
<code>from_xml(element)</code>	Construct a <code>Device</code> object from an XML element.
<code>to_dict([convert_values])</code>	Returns a dictionary containing the contents of the <code>Device</code> object.

`device_type`

Type: `str`

String identifying the type of device.

`classmethod from_dict(d)`

Construct an instance of `Device` from a dictionary.

Parameters `d` (`Mapping[str, Any]`) – The dictionary.

`classmethod from_xml(element)`

Construct a `Device` object from an XML element.

Parameters `element` (`ObjectifiedElement`) – a <Device> XML element from a CEF file

Return type `Device`

`number`

Type: `int`

`to_dict(convert_values=False)`

Returns a dictionary containing the contents of the `Device` object.

Parameters `convert_values` (`bool`) – Recursively convert values into dictionaries, lists etc. as appropriate. Default `False`.

Return type `MutableMapping[str, Any]`

class Flag(*string: str, severity: int*)

Bases: `str`

Represents a flag in a score, to warn that the identification of a compound is poor.

Parameters

- **string** – The text of the flag
- **severity** – The severity of the flag

Methods:

<code>__bool__()</code>	Returns a boolean representation of the <i>Flag</i> .
<code>__eq__(other)</code>	Return <code>self == other</code> .
<code>__ne__(other)</code>	Return <code>self != other</code> .

`__bool__()`

Returns a boolean representation of the *Flag*.

Return type `bool`

`__eq__(other)`

Return `self == other`.

Return type `bool`

`__ne__(other)`

Return `self != other`.

Return type `bool`

typeddict LocationDict

Bases: `TypedDict`

`TypedDict` representing the location of a spectrum within mass spectrometry data.

Optional Keys

- **m** (`float`) – the accurate mass of the compound, determined from the observed mass spectrum.
- **rt** (`float`) – The retention time at which the compound was detected.
- **a** (`float`) – The area of the peak in the EIC.
- **y** (`float`) – The height of the peak in the EIC.

```
class Molecule(name, formula=None, matches=None)
```

Bases: Dictable

Represents a molecule in a CEF file.

Parameters

- **name** (`str`) – The name of the compound
- **formula** (`Union[str, Formula, None]`) – The formula of the compound. If a string it must be parsable by `chemistry_tools.formulae.Formula`. Default `None`.
- **matches** (`Optional[Dict[str, Score]]`) – Dictionary of algo: score match values. Default `None`.

Methods:

<code>__repr__()</code>	Returns a string representation of the <code>Molecule</code> .
<code>__str__()</code>	Returns the molecule as a string.
<code>from_xml(element)</code>	Construct a <code>Molecule</code> object from an XML element.

`__repr__()`

Returns a string representation of the `Molecule`.

Return type `str`

`__str__()`

Returns the molecule as a string.

Return type `str`

`classmethod from_xml(element)`

Construct a `Molecule` object from an XML element.

Parameters `element` (`ObjectifiedElement`) – a Molecule XML element

Return type `Molecule`

`parse_cef(filename)`

Construct an `CompoundList` object from the given .cef file.

Parameters `filename` (`Union[str, Path, PathLike]`) – The filename of the CEF file to read.

Return type `CompoundList`

`parse_compound_scores(element)`

Parse a <CompoundScores> element into a mapping of algorithms to scores.

Parameters `element` (`ObjectifiedElement`) – a CompoundScores XML element.

Return type `Dict[str, Score]`

`parse_match_scores(element)`

Parse a <MatchScores> element into a mapping of algorithms to scores.

Parameters `element` (`ObjectifiedElement`) – a MatchScores XML element.

Return type `Dict[str, Score]`

class `Peak`(*x*, *rx*, *y*, *charge*=0, *label*=")
Bases: `object`

A peak in a Mass Spectrum.

Parameters

- `x` (`float`)
- `rx` (`float`)
- `y` (`float`) – The height of the peak in the EIC.
- `charge` (`int`) – The charge on the peak. Default 0.
- `label` (`str`) – The label of the peak. e.g. “M+H”. Default ''.

Attributes:

<code>charge</code>	The charge on the peak.
<code>label</code>	The label of the peak.
<code>rx</code>	
<code>x</code>	
<code>y</code>	

Methods:

<code>from_dict(d)</code>	Construct an instance of <code>Peak</code> from a dictionary.
<code>from_xml(element)</code>	Construct a <code>Peak</code> object from an XML element.
<code>to_dict([convert_values])</code>	Returns a dictionary containing the contents of the <code>Peak</code> object.

`charge`

Type: `int`

The charge on the peak.

`classmethod from_dict(d)`

Construct an instance of `Peak` from a dictionary.

Parameters `d` (`Mapping[str, Any]`) – The dictionary.

`classmethod from_xml(element)`

Construct a `Peak` object from an XML element.

Parameters `element` (`ObjectifiedElement`) – a <p> XML element from an <MSPeaks> element of a CEF file

Return type `Peak`

`label`

Type: `str`

The label of the peak. e.g. “M+H”

rx**Type:** float**to_dict**(convert_values=False)Returns a dictionary containing the contents of the `Peak` object.**Parameters** convert_values (bool) – Recursively convert values into dictionaries, lists etc. as appropriate. Default `False`.**Return type** MutableMapping[str, Any]**x****Type:** float**y****Type:** float

The height of the peak in the EIC.

class RTRange (start=0.0, end=0.0)Bases: `object`

Represents an <RTRange> element from a CEF file.

Parameters

- **start** (Union[float, timedelta]) – The start time in minutes . Default 0 . 0.
- **end** (Union[float, timedelta]) – The end time in minutes . Default 0 . 0.

Attributes:

<code>end</code>	The end time in minutes
<code>start</code>	The start time in minutes

Methods:

<code>from_dict(d)</code>	Construct an instance of <code>RTRange</code> from a dictionary.
<code>from_xml(element)</code>	Construct ab <code>RTRange</code> object from an XML element.
<code>to_dict([convert_values])</code>	Returns a dictionary containing the contents of the <code>RTRange</code> object.

end**Type:** timedelta

The end time in minutes

classmethod from_dict (d)Construct an instance of `RTRange` from a dictionary.**Parameters** d (Mapping[str, Any]) – The dictionary.**classmethod from_xml (element)**

Construct ab *RTRange* object from an XML element.

Parameters `element` (`ObjectifiedElement`) – The <RTRange> XML element to parse the data from.

Return type `RTRange`

start

Type: `timedelta`

The start time in minutes

to_dict (`convert_values=False`)

Returns a dictionary containing the contents of the *RTRange* object.

Parameters `convert_values` (`bool`) – Recursively convert values into dictionaries, lists etc. as appropriate. Default `False`.

Return type `MutableMapping[str, Any]`

class Score (`score`, `flag_string=''`, `flag_severity=0`)

Bases: `float`

A score indicating how well the compound matches the observed spectrum.

Parameters

- `score` – The score
- `flag_string` (`str`) – Optional flag. See *Flag* for details. Default ''.
- `flag_severity` (`int`) – The severity of the flag. Default 0.

class Spectrum (`spectrum_type=''`, `algorithm=''`, `saturation_limit=0`, `scans=0`, `scan_type=''`,
`ionisation=''`, `polarity=0`, `voltage=0.0`, `device=None`, `peaks=None`, `rt_ranges=None`)

Bases: `Dictable`

Agilent CEF Spectrum.

Parameters

- `spectrum_type` (`str`) – The type of spectrum e.g. 'FbF'. Default ''.
- `algorithm` (`str`) – The algorithm used to identify the compound. Default ''.
- `saturation_limit` (`int`) – Unknown. Might mean saturation limit?. Default 0.
- `scans` (`int`) – Unknown. Presumably the number of scans that make up the spectrum?. Default 0.
- `scan_type` (`str`) – Default ''.
- `ionisation` (`str`) – The type of ionisation e.g. ESI. Default ''.
- `polarity` (`Union[str, int]`) – The polarity of the ionisation. Default 0.
- `device` (`Optional[Device]`) – The device that acquired the data. Default `None`.
- `peaks` (`Optional[Sequence[Peak]]`) – A list of identified peaks in the mass spectrum. Default `None`.
- `rt_ranges` (`Optional[Sequence[RTRange]]`) – A list of retention time ranges for the mass spectrum. Default `None`.

classmethod `from_xml(element)`

Construct a `Spectrum` object from an XML element.

Parameters `element` (`ObjectifiedElement`) – a Spectrum XML element from a CEF file

Return type `Spectrum`

mh_utils.csv_parser

Parser for CSV result files produced by MassHunter Qualitative.

New in version 0.2.0.

<code>ResultParser(raw_results_dir, ...)</code>	Given a directory of CSV results exported from MassHunter, parse them to CSV and JSON.
<code>parse_masshunter_csv(csv_file, csv_outfile, ...)</code>	Parse CSV results files created by MassHunter.

class ResultParser (raw_results_dir, json_results_dir, csv_results_dir)
Bases: `object`

Given a directory of CSV results exported from MassHunter, parse them to CSV and JSON.

Parameters

- **raw_results_dir** (`Union[str, Path, PathLike]`) – The directory in which the raw exports from MassHunter are stored.
- **json_results_dir** (`Union[str, Path, PathLike]`) – The directory to store the output json files in.
- **csv_results_dir** (`Union[str, Path, PathLike]`) – The directory to store the output csv files in.

parse_directory_list (directory_list)

Runs `parse_for_directory()` for each directory in `directory_list`.

Parameters `directory_list` (`Iterable[Union[str, Path, PathLike]]`) – A list of directories to process.

parse_for_directory (directory)

Convert the “CSV Results.csv” file in the given directory to CSV and JSON.

Parameters `directory` (`Union[str, Path, PathLike]`)

parse_masshunter_csv (csv_file, csv_outfile, json_outfile)

Parse CSV results files created by MassHunter.

Parameters

- **csv_file** (`Union[str, Path, PathLike]`)
- **csv_outfile** (`Union[str, Path, PathLike]`)
- **json_outfile** (`Union[str, Path, PathLike]`)

5.3 mh_utils.csv_parser.classes

Classes to model parts of MassHunter CSV files.

New in version 0.2.0.

Classes:

<code>BaseSamplePropertyDict</code>	OrderedDict to store a single property of a set of samples.
<code>Result(cas, name, hits[, index, formula, ...])</code>	Represents a Result in a MassHunter CSV file.
<code>Sample(sample_name, sample_type, ..., [results])</code>	Represents a sample in a MassHunter CSV file.
<code>SampleList([iterable])</code>	A list of <code>mh_utils.csv_parser.classes.Sample</code> objects.
<code>SamplesAreaDict</code>	<code>collections.OrderedDict</code> to store area information parsed from MassHunter results CSV files.
<code>SamplesScoresDict</code>	<code>collections.OrderedDict</code> to store score information parsed from MassHunter results CSV files.

Data:

<code>_R</code>	Invariant TypeVar bound to <code>mh_utils.csv_parser.classes.Result</code> .
<code>_S</code>	Invariant TypeVar bound to <code>mh_utils.csv_parser.classes.Sample</code> .
<code>_SL</code>	Invariant TypeVar bound to <code>mh_utils.csv_parser.classes.SampleList</code> .

class BaseSamplePropertyDict

Bases: `OrderedDict`

OrderedDict to store a single property of a set of samples.

Keys are the sample names and the values are dictionaries mapping compound names to property values.

Attributes:

<code>n_compounds</code>	Returns the number of compounds in the <code>BaseSamplePropertyDict</code> .
<code>n_samples</code>	Returns the number of samples in the <code>BaseSamplePropertyDict</code> .
<code>sample_names</code>	Returns a list of sample names in the <code>BaseSamplePropertyDict</code> .

property n_compounds

Returns the number of compounds in the `BaseSamplePropertyDict`.

Return type `int`

property n_samples

Returns the number of samples in the `BaseSamplePropertyDict`.

Return type `int`

property sample_names

Returns a list of sample names in the `BaseSamplePropertyDict`.

Return type `List[str]`

```
class Result(cas, name, hits, index=-1, formula='', score=0.0, abundance=0, height=0, area=0,
            diff_mDa=0.0, diff_ppm=0.0, rt=0.0, start=0.0, end=0.0, width=0.0, tgt_rt=0.0,
            rt_diff=0.0, mz=0.0, product_mz=0.0, base_peak=0.0, mass=0.0, average_mass=0.0,
            tgt_mass=0.0, mining_algorithm='', z_count=0, max_z=0, min_z=0, n_ions=0,
            polarity='', label='', flags='', flag_severity='', flag_severity_code=0)
```

Bases: `Dictable`

Represents a Result in a MassHunter CSV file.

Parameters

- `cas`
- `name (str)`
- `hits`
- `index (int)` – Default -1.
- `formula (str)` – Default ''.
- `score (float)` – Default 0.0.
- `abundance (float)` – Default 0.
- `height (float)` – Default 0.
- `area (float)` – Default 0.
- `diff_mDa (float)` – Default 0.0.
- `diff_ppm (float)` – Default 0.0.
- `rt (float)` – Default 0.0.
- `start (float)` – Default 0.0.
- `end (float)` – Default 0.0.
- `width (float)` – Default 0.0.
- `tgt_rt (float)` – Default 0.0.
- `rt_diff (float)` – Default 0.0.
- `mz (float)` – Default 0.0.
- `product_mz (float)` – Default 0.0.
- `base_peak (float)` – Default 0.0.
- `mass (float)` – Default 0.0.
- `average_mass (float)` – Default 0.0.
- `tgt_mass (float)` – Default 0.0.
- `mining_algorithm (str)` – Default ''.
- `z_count (int)` – Default 0.
- `max_z (int)` – Default 0.
- `min_z (int)` – Default 0.
- `n_ions (int)` – Default 0.
- `polarity (str)` – Default ''.
- `label (str)` – Default ''.
- `flags (str)` – Default ''.
- `flag_severity (str)` – Default ''.
- `flag_severity_code (int)` – Default 0.

Methods:

<code>from_series(series)</code>	Construct a <code>Result</code> from a <code>pandas.Series</code> .
----------------------------------	---

classmethod from_series(series)

Construct a `Result` from a `pandas.Series`.

Parameters `series (Series)`

Return type `~R`

```
class Sample(sample_name, sample_type, instrument_name, position, user, acq_method, da_method,
             irm_cal_status, filename, results=None)
```

Bases: Dictable

Represents a sample in a MassHunter CSV file.

Parameters

- `sample_name`
- `sample_type`
- `instrument_name`
- `position`
- `user`
- `acq_method`
- `da_method`
- `irm_cal_status`
- `filename`
- `results` – Default `None`.

Methods:

<code>add_result(result)</code>	Add a result to the sample.
<code>from_series(series)</code>	Construct a <code>Sample</code> from a <code>pandas.Series</code> .

Attributes:

<code>results_list</code>	Returns a list of results in the order in which they were identified.
---------------------------	---

`add_result(result)`

Add a result to the sample.

Parameters `result`

`classmethod from_series(series)`

Construct a `Sample` from a `pandas.Series`.

Parameters `series`

Return type `~S`

Returns

`property results_list`

Returns a list of results in the order in which they were identified.

I.e. sorted by the `Cpd` value from the csv export.

Return type `List[Result]`

class SampleList (iterable=(), /)
Bases: `List[Sample]`
A list of `mh_utils.csv_parser.classes.Sample` objects.

Methods:

<code>add_new_sample(*args, **kwargs)</code>	Add a new sample to the list and return the <code>Sample</code> object representing it.
<code>add_sample(sample)</code>	Add a <code>Sample</code> object to the list.
<code>add_sample_from_series(series)</code>	Create a new sample object from a <code>pandas.Series</code> and add it to the list.
<code>filter(sample_names[, key, exclude])</code>	Filter the list to only contain <code>sample_names</code> whose name is in <code>sample_names</code> .
<code>from_json_file(filename, **kwargs)</code>	Construct a <code>SampleList</code> from JSON file.
<code>get_areas_and_scores(compound_name[, ...])</code>	Returns two dictionaries: one containing sample names and peak areas for the compound with the given name, the other containing sample names and scores.
<code>get_areas_and_scores_for_compounds(...)</code>	Returns two dictionaries: one containing sample names and peak areas for the compounds with the given names, the other containing sample names and scores.
<code>get_areas_for_compounds(compound_names[, ...])</code>	Returns a dictionary containing sample names and peak areas for the compounds with the given names.
<code>get_compounds()</code>	Returns a list containing the names of the compounds present in the samples in alphabetical order.
<code>get_peak_areas(compound_name[, include_none])</code>	Returns a dictionary containing sample names and peak areas for the compound with the given name.
<code>get_retention_times(compound_name[, ...])</code>	Returns a dictionary containing sample names and retention times for the compound with the given name.
<code>get_scores(compound_name[, include_none])</code>	Returns a dictionary containing sample names and scores for the compound with the given name.
<code>rename_samples(rename_mapping[, key])</code>	Rename the samples in the list.
<code>reorder_samples(order_mapping[, key])</code>	Reorder the list of <code>Samples</code> in place.
<code>sort_samples(key[, reverse])</code>	Sort the list of <code>Samples</code> in place.

Attributes:

<code>sample_names</code>	Returns a list of sample names in the <code>SampleList</code> .
---------------------------	---

`add_new_sample(*args, **kwargs)`
Add a new sample to the list and return the `Sample` object representing it.

`add_sample(sample)`
Add a `Sample` object to the list.

Parameters `sample (Sample)`

Return type `Sample`

add_sample_from_series (*series*)

Create a new sample object from a pandas.Series and add it to the list.

Return type *Sample*

Returns The newly created Sample object.

Parameters *series* (Series)

filter (*sample_names*, *key='sample_name'*, *exclude=False*)

Filter the list to only contain sample_names whose name is in sample_names.

Parameters

- **sample_names** (Iterable[str]) – A list of sample names to include
- **key** (str) – The name of the property in the sample to sort by. Default 'sample_name'.
- **exclude** (bool) – If True, any sample whose name is in sample_names will be excluded from the output, rather than included. Default False.

Return type *~SL*

classmethod from_json_file (*filename*, ***kwargs*)

Construct a *SampleList* from JSON file.

Parameters

- **filename** (Union[str, Path, PathLike]) – The filename of the JSON file.
- ****kwargs** – Keyword arguments passed to `domdf_python_tools.paths.PathPlus.load_json()`.

Return type *~SL*

get_areas_and_scores (*compound_name*, *include_none=False*)

Returns two dictionaries: one containing sample names and peak areas for the compound with the given name, the other containing sample names and scores.

Parameters

- **compound_name** (str)
- **include_none** (bool) – Whether samples where the compound was not found should be included in the results. Default False.

Return type Tuple[OrderedDict, OrderedDict]

get_areas_and_scores_for_compounds (*compound_names*, *include_none=False*)

Returns two dictionaries: one containing sample names and peak areas for the compounds with the given names, the other containing sample names and scores.

Parameters

- **compound_names** (Iterable[str])
- **include_none** (bool) – Whether samples where none of the specified compounds were found should be included in the results. Default False.

Return type Tuple[SamplesAreaDict, SamplesScoresDict]

get_areas_for_compounds (*compound_names*, *include_none=False*)

Returns a dictionary containing sample names and peak areas for the compounds with the given names.

Parameters

- **compound_names** (`Iterable[str]`)
- **include_none** (`bool`) – Whether samples where none of the specified compounds were found should be included in the results. Default `False`.

Return type `SamplesAreaDict`**get_compounds()**

Returns a list containing the names of the compounds present in the samples in alphabetical order.

Return type `List[str]`**get_peak_areas** (*compound_name*, *include_none=False*)

Returns a dictionary containing sample names and peak areas for the compound with the given name.

Parameters

- **compound_name** (`str`)
- **include_none** (`bool`) – Whether samples where the compound was not found should be included in the results. Default `False`.

Return type `OrderedDict`**get_retention_times** (*compound_name*, *include_none=False*)

Returns a dictionary containing sample names and retention times for the compound with the given name.

Parameters

- **compound_name** (`str`)
- **include_none** (`bool`) – Whether samples where the compound was not found should be included in the results. Default `False`.

Return type `OrderedDict`**get_scores** (*compound_name*, *include_none=False*)

Returns a dictionary containing sample names and scores for the compound with the given name.

Parameters

- **compound_name** (`str`)
- **include_none** (`bool`) – Whether samples where the compound was not found should be included in the results. Default `False`.

Return type `OrderedDict`

rename_samples (*rename_mapping*, *key='sample_name'*)

Rename the samples in the list.

Parameters

- **rename_mapping** (`Dict`) – A mapping between current sample names and their new names.
- **key** (`str`) – The name of the property in the sample to sort by. Default '`sample_name`'.

Use `rename_mapping=:`py:`obj=None` or omit the sample from the `rename_mapping` entirely to leave the name unchanged.

For example:

```
rename_mapping = {
    "Propellant lug +ve": "Alliant Unique 1µg/L +ESI",
    "Propellant 1mg +ve": "Alliant Unique 1mg/L +ESI",
    "Propellant 1mg -ve": None,
}
```

reorder_samples (*order_mapping*, *key='sample_name'*)

Reorder the list of Samples in place.

Parameters

- **order_mapping** (`Dict`) – A mapping between sample names and their new position in the list. For example:

```
order_mapping = {
    "Propellant lug +ve": 0,
    "Propellant 1mg +ve": 1,
    "Propellant lug -ve": 2,
    "Propellant 1mg -ve": 3,
}
```

- **key** (`str`) – The name of the property in the sample to sort by. Default '`sample_name`'.

property sample_names

Returns a list of sample names in the *SampleList*.

Return type `List[str]`

sort_samples (*key*, *reverse=False*)

Sort the list of Samples in place.

Parameters

- **key** (`str`) – The name of the property in the sample to sort by.
- **reverse** (`bool`) – Whether the list should be sorted in reverse order. Default `False`.

Return type

class SamplesAreaDictBases: *BaseSamplePropertyDict*

collections.OrderedDict to store area information parsed from MassHunter results CSV files.

Methods:

<code>get_compound_areas(compound_name)</code>	Get the peak areas for the given compound in every sample.
--	--

get_compound_areas (compound_name)

Get the peak areas for the given compound in every sample.

Parameters `compound_name` (`str`)**Return type** `List[float]`**class SamplesScoresDict**Bases: *BaseSamplePropertyDict*

collections.OrderedDict to store score information parsed from MassHunter results CSV files.

Methods:

<code>get_compound_scores(compound_name)</code>	Get the peak scores for the given compound in every sample.
---	---

get_compound_scores (compound_name)

Get the peak scores for the given compound in every sample.

Parameters `compound_name` (`str`)**Return type** `List[float]`**_R = TypeVar(_R, bound=Result)****Type:** `TypeVar`Invariant `TypeVar` bound to `mh_utils.csv_parser.classes.Result`.**_S = TypeVar(_S, bound=Sample)****Type:** `TypeVar`Invariant `TypeVar` bound to `mh_utils.csv_parser.classes.Sample`.**_SL = TypeVar(_SL, bound=SampleList)****Type:** `TypeVar`Invariant `TypeVar` bound to `mh_utils.csv_parser.classes.SampleList`.

5.4 mh_utils.csv_parser.utils

CSV utility functions.

New in version 0.2.0.

Functions:

<code>concatenate_json(*files[, outfile])</code>	Concatenate multiple JSON files together and return a list of Sample objects in the concatenated json output.
<code>drop_columns(df, *[, axis, inplace])</code>	Drop columns from the MassHunter CSV file.
<code>reorder_columns(df)</code>	Reorder columns from the MassHunter CSV file.

`concatenate_json(*files, outfile=None)`

Concatenate multiple JSON files together and return a list of Sample objects in the concatenated json output.

Parameters

- `*files` (`Union[str, Path, PathLike]`) – The files to concatenate.
- `outfile` (`Union[str, Path, PathLike, None]`) – The file to save the output as. If `None` no file will be saved. Default `None`.

Return type `SampleList`

`drop_columns(df, *, axis=1, inplace=True, **kwargs)`

Drop columns from the MassHunter CSV file.

Parameters

- `df` (`DataFrame`) – The `pandas.DataFrame` to drop columns in.
- `axis` (`int`) – Which axis to drop columns on. Default 1.
- `inplace` (`bool`) – Whether to modify the `pandas.DataFrame` in place. Default `True`.
- `kwargs` – Additional keyword arguments passed to `pandas.DataFrame.drop()`.

Return type `DataFrame`

`reorder_columns(df)`

Reorder columns from the MassHunter CSV file.

Parameters `df` (`DataFrame`) – The `pandas.DataFrame` to reorder columns in.

Return type `DataFrame`

mh_utils.worklist_parser

Parser for MassHunter worklists.

Only one function is defined here: `read_worklist()`, which reads the given worklist file and returns a `mh_utils.worklist_parser.classes.Worklist` file representing it. The other functions and classes must be imported from submodules of this package.

read_worklist(filename)

Read the worklist from the given file.

Parameters `filename` (`Union[str, Path, PathLike]`) – The filename of the worklist.

Return type `Worklist`

6.2 mh_utils.worklist_parser.classes

Main classes for the worklist parser.

Classes:

<code>Attribute(attribute_id, attribute_type, ...)</code>	Represents an Attribute.
<code>Checksum(SchemaVersion, ALGO_VERSION, HASHCODE)</code>	Represents a checksum for a worklist.
<code>JobData(id, job_type, run_status[, sample_info])</code>	Represents an entry in the worklist.
<code>Macro(project_name, procedure_name, ...)</code>	Represents a macro in a worklist.
<code>Worklist(version, locked_run_mode, ...)</code>	Class that represents an Agilent MassHunter worklist.

class Attribute(attribute_id, attribute_type, field_type, system_name, header_name, data_type, default_data_value, reorder_id, show_hide_status, column_width)

Represents an Attribute.

Parameters

- **attribute_id** (`int`)
- **attribute_type** (`AttributeType`)
 - The attribute type identifier.
- **field_type** (`int`) – The field type identifier.
- **system_name** (`str`)
- **header_name** (`str`)
- **data_type** (`int`)
- **default_data_value** (`str`)
- **reorder_id** (`int`)
- **show_hide_status** (`Union[str, bool]`)
- **column_width** (`int`) – end{multicols}

Attributes:

<code>attribute_id</code>	
<code>attribute_type</code>	The attribute type identifier.
<code>column_width</code>	
<code>data_type</code>	
<code>default_data_value</code>	
<code>field_type</code>	The field type identifier.
<code>header_name</code>	
<code>reorder_id</code>	
<code>show_hide_status</code>	
<code>system_name</code>	

Methods:

<code>from_dict(d)</code>	Construct an instance of <code>Attribute</code> from a dictionary.
<code>from_xml(element)</code>	Construct an <code>Attribute</code> object from an XML element.
<code>to_dict([convert_values])</code>	Returns a dictionary containing the contents of the <code>Attribute</code> object.

attribute_id
Type: `int`

attribute_type
Type: `AttributeType`

The attribute type identifier.

Can be System Defined (0), System Used (1), or User Added (2).

column_width
Type: `int`

data_type
Type: `Any`

default_data_value
Type: `str`

field_type
Type: `int`

The field type identifier.

Each of the system defined columns have a field type starting from sampleid = 0 to reserved6 = 24.

The system used column can be ‘compound param’ = 35, ‘optim param’ = 36, ‘mass param’ = 37 and ‘protein param’ = 38.

The User added columns start from 45.

classmethod from_dict(d)

Construct an instance of `Attribute` from a dictionary.

Parameters `d` (`Mapping[str, Any]`) – The dictionary.

classmethod from_xml(element)

Construct an `Attribute` object from an XML element.

Return type `Attribute`

header_name

Type: `str`

reorder_id

Type: `int`

show_hide_status

Type: `bool`

system_name

Type: `str`

to_dict(convert_values=False)

Returns a dictionary containing the contents of the `Attribute` object.

Parameters `convert_values` (`bool`) –

Recursively convert values into dictionaries, lists etc. as appropriate. Default `False`.

Return type `MutableMapping[str, Any]`

class Checksum(SchemaVersion, ALGO_VERSION, HASHCODE)

Represents a checksum for a worklist.

The format of the checksum is unknown.

Parameters

- `SchemaVersion` (`int`)
- `ALGO_VERSION` (`int`)
- `HASHCODE` (`str`)

Attributes:

`ALGO_VERSION`

`HASHCODE`

`SchemaVersion`

Methods:

<code>from_dict(d)</code>	Construct an instance of <code>Checksum</code> from a dictionary.
<code>from_xml(element)</code>	Construct a <code>Checksum</code> object from an XML element.
<code>to_dict([convert_values])</code>	Returns a dictionary containing the contents of the <code>Checksum</code> object.

ALGO_VERSION**Type:** int**HASHCODE****Type:** str**SchemaVersion****Type:** int**classmethod from_dict(d)**Construct an instance of `Checksum` from a dictionary.**Parameters** `d` (Mapping[str, Any]) – The dictionary.**classmethod from_xml(element)**Construct a `Checksum` object from an XML element.**Return type** `Checksum`**to_dict(convert_values=False)**Returns a dictionary containing the contents of the `Checksum` object.**Parameters** `convert_values` (bool) – Recursively convert values into dictionaries, lists etc. as appropriate. Default `False`.**Return type** MutableMapping[str, Any]**class JobData(id, job_type, run_status, sample_info=None)****Bases:** Dictable

Represents an entry in the worklist.

Parameters

- `id` (Union[str, UUID]) – The ID of the job.
- `job_type` (int) – The type of job. TODO: enum of values
- `run_status` (int) – The status of the analysis. TODO: enum of values
- `sample_info` (Optional[dict]) – Optional key: value mapping of information about the sample. Default `None`.

Methods:

<code>from_xml(element[, user_columns])</code>	Construct a <code>JobData</code> object from an XML element.
--	--

classmethod from_xml(element, user_columns=None)

Construct a `JobData` object from an XML element.

Parameters

- **element** (`ObjectifiedElement`) – The XML element to parse the data from
- **user_columns** (`Optional[Dict[str, Column]]`) – Optional mapping of user column labels to `Column` objects. Default `None`.

Return type `JobData`

```
class Macro(project_name, procedure_name, input_parameter, output_data_type,
               output_parameter, display_string)
```

Represents a macro in a worklist.

Parameters

- **project_name** (`str`)
- **procedure_name** (`str`)
- **input_parameter** (`str`)
- **output_data_type** (`int`)
- **output_parameter** (`str`) – .
- **display_string** (`str`)

Attributes:

<i>display_string</i>	
<i>input_parameter</i>	
<i>output_data_type</i>	
<i>output_parameter</i>	
<i>procedure_name</i>	
<i>project_name</i>	
<i>undefined</i>	Returns whether the macro is undefined.

Methods:

<code>from_dict(d)</code>	Construct an instance of <code>Macro</code> from a dictionary.
<code>from_xml(element)</code>	Construct a <code>Macro</code> object from an XML element.
<code>to_dict([convert_values])</code>	Returns a dictionary containing the contents of the <code>Macro</code> object.

display_string
Type: `str`

classmethod from_dict(d)
Construct an instance of `Macro` from a dictionary.

Parameters `d` (`Mapping[str, Any]`) – The dictionary.

classmethod from_xml(element)
Construct a `Macro` object from an XML element.

Return type `Macro`

input_parameter
Type: `str`

output_data_type
Type: `int`

output_parameter

Type: str

procedure_name

Type: str

project_name

Type: str

to_dict (convert_values=False)Returns a dictionary containing the contents of the *Macro* object.**Parameters** convert_values (bool) – Recursively convert values into dictionaries, lists etc. as appropriate. Default False.**Return type** MutableMapping[str, Any]**property undefined**

Returns whether the macro is undefined.

Return type bool**class Worklist** (version, locked_run_mode, instrument_name, params, user_columns, jobs, checksum)

Bases: XMLFileMixin, Dictable

Class that represents an Agilent MassHunter worklist.

Parameters

- **version** (float) – WorklistInfo version number
- **locked_run_mode** (bool) – Flag to indicate whether the data was acquired in locked mode. Yes = -1. No = 0.
- **instrument_name** (str) – The name of the instrument.
- **params** (dict) – Mapping of parameter names to values. TODO: Check
- **user_columns** (dict) – Mapping of user columns to ??? TODO
- **jobs** (Sequence[JobData])
- **checksum** (Checksum) – The checksum of the worklist file. The format is unknown.

Methods:

__repr__()	Return a string representation of the <i>Worklist</i> .
as_dataframe()	Returns the <i>Worklist</i> as a pandas. DataFrame.
from_xml(element)	Construct a <i>Worklist</i> object from an XML element.

__repr__()Return a string representation of the *Worklist*.**Return type** str**as_dataframe()**Returns the *Worklist* as a pandas. DataFrame.

Return type `DataFrame`

classmethod from_xml(element)
Construct a `Worklist` object from an XML element.

Return type `Worklist`

6.3 mh_utils.worklist_parser.columns

Properties for columns in a Worklist.

<code>Column(name, attribute_id, attribute_type, ...)</code>	Represents a column in a worklist.
<code>injection_volume(val)</code>	Handle special case for injection volume of -1, which indicates “As Method”.
<code>columns</code>	Mapping of column names to column objects.

class Column(name, attribute_id, attribute_type, dtype, default_value, field_type=None, reorder_id=None)

Bases: `object`

Represents a column in a worklist.

Parameters

- **name** (`str`) – The name of the column
- **attribute_id** (`int`)
- **attribute_type** (`AttributeType`) – The attribute type identifier.
- **dtype** (`Callable`) – The field datatype.
- **default_value** (`Any`)
- **field_type** (`Optional[int]`) – The field type identifier. Default `None`.
- **reorder_id** (`Optional[int]`) – Default `None`.

attribute_id

Type: `int`

attribute_type

Type: `AttributeType`

The attribute type identifier.

Can be System Defined (0), System Used (1), or User Added (2).

cast_value(value)

Cast value to the dtype of this column.

default_value

Type: `Any`

dtype

Type: `Callable`

field_type

Type: `Optional[int]`

The field type identifier.

Each of the system defined columns have a field type starting from sampleid = 0 to reserved6 = 24.

The system used column can be ‘compound param’ = 35, ‘optim param’ = 36, ‘mass param’ = 37 and ‘protein param’ = 38.

The User added columns start from 45.

classmethod from_attribute(*attribute*)

Construct a column for a `mh_utils.worklist_parser.classes.Attribute`.

Parameters `attribute (Attribute) – Attribute`

Return type `Column`

classmethod from_dict(*d*)

Construct an instance of `Column` from a dictionary.

Parameters `d (Mapping[str, Any]) – The dictionary.`

name

Type: `str`

The name of the column

reorder_id

Type: `Optional[int]`

to_dict(*convert_values=False*)

Returns a dictionary containing the contents of the `Column` object.

Parameters `convert_values (bool) – Recursively convert values into dictionaries, lists etc. as appropriate. Default False.`

Return type `MutableMapping[str, Any]`

injection_volume(*val*)

Handle special case for injection volume of -1, which indicates “As Method”.

Parameters `val (Union[float, str])`

Returns

Return type `Union[int, str]`

columns

Mapping of column names to column objects.

6.4 mh_utils.worklist_parser.enums

Enumerations of values.

Classes:

<code>AttributeType(value)</code>	Enumeration of values for column/attribute types.
-----------------------------------	---

enum AttributeType (value)

Bases: `enum_tools.custom_enums.IntEnum`

Enumeration of values for column/attribute types.

Member Type `int`

Valid values are as follows:

`SystemDefined = <AttributeType.SystemDefined: 0>`

`SystemUsed = <AttributeType.SystemUsed: 1>`

`UserAdded = <AttributeType.UserAdded: 2>`

6.5 mh_utils.worklist_parser.parser

MassHunter worklist parser.

Functions:

<code>parse_datetime(the_date)</code>	Parse a datetime from a worklist or contents file.
<code>parse_params(element)</code>	Parse the worklist execution parameters from XML.
<code>parse_sample_info(element[, user_columns])</code>	Parse information about a sample in a worklist from XML.

Data:

<code>sample_info_tags</code>	Mapping of XML tag names to attribute names.
-------------------------------	--

parse_datetime (the_date)

Parse a datetime from a worklist or contents file.

Parameters `the_date (str)` – The date and time as a string in the format

`%Y-%m-%dT%H:%M:%S%z`.

Return type `datetime`

parse_params (element)

Parse the worklist execution parameters from XML.

Parameters `element (ObjectifiedElement)`

Return type `Dict[str, Any]`

Returns Mapping of keys to parameter values.

parse_sample_info (*element*, *user_columns=None*)
Parse information about a sample in a worklist from XML.

Parameters

- **element** (*ObjectifiedElement*) – The XML element to parse the data from
- **user_columns** (*Optional[Dict[str, Column]]*) – Optional mapping of user column labels to *Column* objects. Default *None*.

Return type *Dict[str, Any]*

sample_info_tags
Mapping of XML tag names to attribute names.

6.6 Example Usage

Listing 1: *read_worklist.py*

```

1 # 3rd party
2 from pandas import DataFrame
3
4 # this package
5 from mh_utils.worklist_parser import Worklist, read_worklist
6
7 # Replace 'worklist.wkl' with the filename of your worklist.
8 wkl: Worklist = read_worklist("worklist.wkl")
9
10 df: DataFrame = wkl.as_dataframe()
11
12 # Filter columns
13 df = df[["Sample Name", "Method", "Data File", "Drying Gas", "Gas Temp",
14          ↪"Nebulizer"]]
15
16 # Get just the filename from 'Method' and 'Data File'
17 df["Method"] = [x.name for x in df["Method"]]
18 df["Data File"] = [x.name for x in df["Data File"]]
19
20 # Show the DataFrame
21 print(df.to_string())

```

Listing 2: *worklist_df_head.txt*

	Sample Name	Method	
0	Methanol Blank +ve	Maitre Gunshot Residue Positive.m	↳
	↳Methanol_Blank_+ve_191121-0001-r001.d		↳
1	Propellant 1mg +ve	Maitre Gunshot Residue Positive.m	↳
	↳Propellant_1mg_+ve_191121-0002-r001.d		↳
2	Propellant lug +ve	Maitre Gunshot Residue Positive.m	↳
	↳Propellant_lug_+ve_191121-0003-r001.d		↳
3	Methanol Blank -ve	Maitre Gunshot Residue Negative.m	↳
	↳Methanol_Blank_-ve_191121-0004-r001.d		↳

(continues on next page)

(continued from previous page)

4	Propellant 1mg -ve	Maitre Gunshot Residue Negative.m	✓
↳	Propellant_1mg_-ve_191121-0005-r001.d		
5	Propellant 1ug -ve	Maitre Gunshot Residue Negative.m	✓
↳	Propellant_1ug_-ve_191121-0006-r001.d		
6	Methanol Blank +ve 5ul	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Methanol_Blank_+ve_5ul_191121-0007-r001.d		
7	Propellant 1mg +ve 5ul	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_+ve_5ul_191121-0008-r001.d		
8	Propellant 1ug +ve 5ul	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_+ve_5ul_191121-0009-r001.d		
9	Methanol Blank	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Methanol_Blank_191121-0010-r001.d		
10	Propellant 1mg gas 200	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_gas_200_191121-0011-r001.d	200	
11	Propellant 1ug gas 200	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_gas_200_191121-0012-r001.d	200	
12	Propellant 1mg gas 280	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_gas_280_191121-0013-r001.d	280	
13	Propellant 1ug gas 280	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_gas_280_191121-0014-r001.d	280	
14	Propellant 1mg drying 14	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_drying_14_191121-0015-r001.d	14	
15	Propellant 1ug drying 14	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_drying_14_191121-0016-r001.d	14	
16	Propellant 1mg drying 16	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_drying_16_191121-0017-r001.d	16	
17	Propellant 1ug drying 16	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_drying_16_191121-0018-r001.d	16	
18	Propellant 1mg drying 18	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_drying_18_191121-0019-r001.d	18	
19	Propellant 1ug drying 18	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_drying_18_191121-0020-r001.d	18	
20	Propellant 1mg nebul 40	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_nebul_40_191121-0021-r001.d	40	
21	Propellant 1ug nebul 40	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_nebul_40_191121-0022-r001.d	40	
22	Propellant 1mg nebul 50	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_nebul_50_191121-0023-r001.d	50	
23	Propellant 1ug nebul 50	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_nebul_50_191121-0024-r001.d	50	
24	Propellant 1mg nebul 60	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1mg_nebul_60_191121-0025-r001.d	60	
25	Propellant 1ug nebul 60	Maitre Gunshot Residue Positive 5ul.m	✓
↳	Propellant_1ug_nebul_60_191121-0026-r001.d	60	
...			

Listing 3: read_worklist.py

```
21 # save as CSV
22 df.to_csv("worklist.csv")
24
25 # save as JSON
26 df.to_json("worklist.json", indent=2)
```

6.6.1 Output

Python Module Index

m

`mh_utils.cef_parser`, 9
`mh_utils.csv_parser`, 21
`mh_utils.csv_parser.classes`, 22
`mh_utils.csv_parser.utils`, 30
`mh_utils.utils`, 5
`mh_utils.worklist_parser`, 31
`mh_utils.worklist_parser.columns`, 41
`mh_utils.worklist_parser.enums`, 43
`mh_utils.worklist_parser.parser`, 43
`mh_utils.xml`, 7

Index

Symbols

_R (*in module mh_utils.csv_parser.classes*), 29
_S (*in module mh_utils.csv_parser.classes*), 29
_SL (*in module mh_utils.csv_parser.classes*), 29
__bool__() (*Flag method*), 14
__eq__() (*Flag method*), 14
__ne__() (*Flag method*), 14
__repr__() (*Compound method*), 11
__repr__() (*CompoundList method*), 12
__repr__() (*Molecule method*), 15
__repr__() (*Worklist method*), 39
__str__() (*Compound method*), 11
__str__() (*CompoundList method*), 12
__str__() (*Molecule method*), 15

A

add_new_sample() (*SampleList method*), 25
add_result() (*Sample method*), 24
add_sample() (*SampleList method*), 25
add_sample_from_series() (*SampleList method*), 26
algo (*Compound attribute*), 11
ALGO_VERSION (*Checksum attribute*), 35
as_dataframe() (*Worklist method*), 39
as_path() (*in module mh_utils.utils*), 5
Attribute (*class in mh_utils.worklist_parser.classes*), 31
attribute_id (*Attribute attribute*), 32
attribute_id (*Column attribute*), 41
attribute_type (*Attribute attribute*), 32
attribute_type (*Column attribute*), 41

B

BaseSamplePropertyDict (*class in mh_utils.csv_parser.classes*), 22

C

camel_to_snake() (*in module mh_utils.utils*), 5
cast_value() (*Column method*), 41
charge (*Peak attribute*), 16
Checksum (*class in mh_utils.worklist_parser.classes*), 34
Column (*class in mh_utils.worklist_parser.columns*), 41
column_width (*Attribute attribute*), 32
columns (*in module mh_utils.worklist_parser.columns*), 42
Compound (*class in mh_utils.cef_parser*), 10
compound_scores (*Compound attribute*), 11

CompoundList (*class in mh_utils.cef_parser*), 12
concatenate_json() (*in module mh_utils.csv_parser.utils*), 30

D

data_type (*Attribute attribute*), 32
default_data_value (*Attribute attribute*), 32
default_value (*Column attribute*), 41
Device (*class in mh_utils.cef_parser*), 13
device_type (*Device attribute*), 13
display_string (*Macro attribute*), 37
drop_columns() (*in module mh_utils.csv_parser.utils*), 30
dtype (*Column attribute*), 41

E

element_to_bool() (*in module mh_utils.utils*), 5
end (*RTRange attribute*), 17

F

field_type (*Attribute attribute*), 32
field_type (*Column attribute*), 42
filter() (*SampleList method*), 26
Flag (*class in mh_utils.cef_parser*), 14
from_attribute() (*Column class method*), 42
from_dict() (*Attribute class method*), 34
from_dict() (*Checksum class method*), 35
from_dict() (*Column class method*), 42
from_dict() (*Device class method*), 13
from_dict() (*Macro class method*), 37
from_dict() (*Peak class method*), 16
from_dict() (*RTRange class method*), 17
from_json_file() (*SampleList class method*), 26
from_series() (*Result class method*), 23
from_series() (*Sample class method*), 24
from_xml() (*Attribute class method*), 34
from_xml() (*Checksum class method*), 35
from_xml() (*Compound class method*), 11
from_xml() (*CompoundList class method*), 12
from_xml() (*Device class method*), 13
from_xml() (*JobData class method*), 35
from_xml() (*Macro class method*), 37
from_xml() (*Molecule class method*), 15
from_xml() (*Peak class method*), 16
from_xml() (*RTRange class method*), 17
from_xml() (*Spectrum class method*), 18

from_xml () (*Worklist class method*), 41
from_xml () (*XMLFileMixin class method*), 7
from_xml_file () (*XMLFileMixin class method*), 7

G

get_areas_and_scores () (*SampleList method*), 26
get_areas_and_scores_for_compounds () (*SampleList method*), 26
get_areas_for_compounds () (*SampleList method*), 27
get_compound_areas () (*SamplesAreaDict method*), 29
get_compound_scores () (*SamplesScoresDict method*), 29
get_compounds () (*SampleList method*), 27
get_peak_areas () (*SampleList method*), 27
get_retention_times () (*SampleList method*), 27
get_scores () (*SampleList method*), 27
get_validated_tree () (*in module mh_utils.xml*), 7

module, 43
mh_utils.xml
 module, 7
module
 mh_utils.cef_parser, 9
 mh_utils.csv_parser, 21
mh_utils.csv_parser.classes, 22
mh_utils.csv_parser.utils, 30
mh_utils.utils, 5
mh_utils.worklist_parser, 31
 mh_utils.worklist_parser.columns, 41
 mh_utils.worklist_parser.enums, 43
 mh_utils.worklist_parser.parser, 43
 mh_utils.xml, 7
Molecule (*class in mh_utils.cef_parser*), 15

H

HASHCODE (*Checksum attribute*), 35
header_name (*Attribute attribute*), 34

I

injection_volume () (*in module mh_utils.worklist_parser.columns*), 42
input_parameter (*Macro attribute*), 37
instrument (*CompoundList attribute*), 12

J

JobData (*class in mh_utils.worklist_parser.classes*), 35

L

label (*Peak attribute*), 16
location (*Compound attribute*), 11
LocationDict (*typeddict in mh_utils.cef_parser*), 14

M

Macro (*class in mh_utils.worklist_parser.classes*), 37
make_timedelta () (*in module mh_utils.utils*), 5
mh_utils.cef_parser
 module, 9
mh_utils.csv_parser
 module, 21
mh_utils.csv_parser.classes
 module, 22
mh_utils.csv_parser.utils
 module, 30
mh_utils.utils
 module, 5
mh_utils.worklist_parser
 module, 31
mh_utils.worklist_parser.columns
 module, 41
mh_utils.worklist_parser.enums
 module, 43
mh_utils.worklist_parser.parser

n_compounds () (*BaseSamplePropertyDict property*), 22
n_samples () (*BaseSamplePropertyDict property*), 22
name (*Column attribute*), 42
number (*Device attribute*), 13

output_data_type (*Macro attribute*), 37
output_parameter (*Macro attribute*), 37

P

parse_cef () (*in module mh_utils.cef_parser*), 15
parse_compound_scores () (*in module mh_utils.cef_parser*), 1
parse_datetime () (*in module mh_utils.worklist_parser.parser*),
parse_directory_list () (*ResultParser method*), 21
parse_for_directory () (*ResultParser method*), 21
parse_masshunter_csv () (*in module mh_utils.csv_parser*), 21
parse_match_scores () (*in module mh_utils.cef_parser*), 15
parse_params () (*in module mh_utils.worklist_parser.parser*), 43
parse_sample_info () (*in module mh_utils.worklist_parser.parser*), 21
Peak (*class in mh_utils.cef_parser*), 16
procedure_name (*Macro attribute*), 39
project_name (*Macro attribute*), 39

R

read_worklist () (*in module mh_utils.worklist_parser*), 31
rename_samples () (*SampleList method*), 28
reorder_columns () (*in module mh_utils.csv_parser.utils*), 30
reorder_id (*Attribute attribute*), 34
reorder_id (*Column attribute*), 42
reorder_samples () (*SampleList method*), 28
Result (*class in mh_utils.csv_parser.classes*), 23
ResultParser (*class in mh_utils.csv_parser*), 21
results (*Compound attribute*), 11
results_list () (*Sample property*), 24
RTRange (*class in mh_utils.cef_parser*), 17
rx (*Peak attribute*), 16

S

Sample (*class in mh_utils.csv_parser.classes*), 24
sample_info_tags (*in module mh_utils.worklist_parser.parser*), 44
sample_names () (*BaseSamplePropertyDict property*), 23
sample_names () (*SampleList property*), 28
SampleList (*class in mh_utils.csv_parser.classes*), 25
SamplesAreaDict (*class in mh_utils.csv_parser.classes*), 29
SamplesScoresDict (*class in mh_utils.csv_parser.classes*), 29
SchemaVersion (*Checksum attribute*), 35
Score (*class in mh_utils.cef_parser*), 18
show_hide_status (*Attribute attribute*), 34
sort_samples () (*SampleList method*), 28
spectra (*Compound attribute*), 11
Spectrum (*class in mh_utils.cef_parser*), 18
start (*RTRange attribute*), 18
strip_string () (*in module mh_utils.utils*), 6
system_name (*Attribute attribute*), 34
SystemDefined (*AttributeType attribute*), 43
SystemUsed (*AttributeType attribute*), 43

T

to_dict () (*Attribute method*), 34
to_dict () (*Checksum method*), 35
to_dict () (*Column method*), 42
to_dict () (*Device method*), 13
to_dict () (*Macro method*), 39
to_dict () (*Peak method*), 17
to_dict () (*RTRange method*), 18

U

undefined () (*Macro property*), 39
UserAdded (*AttributeType attribute*), 43

W

Worklist (*class in mh_utils.worklist_parser.classes*), 39

X

x (*Peak attribute*), 17
XMLFileMixin (*class in mh_utils.xml*), 7

Y

y (*Peak attribute*), 17